
Software Release Note

IMS Wf & SC Ultra

Release 3.7 Patch #10

Release Date

Date	Release	Applicability
September 3, 2025	3.7.10	This software release applies to <ul style="list-style-type: none">Correction R1871

Modifications in Software Release 3.7.10 for IMS Wf & SC Ultra

1. List of Bugs Fixed

ID	Application	Description
1871	Acquisition	Fixes a bug with lost CKB data. When the bug was present, when a user restarted a CheckerBoard acquisition and clicked Cancel in the Save As window, the data from the previous acquisition was lost. This has been corrected.